



Product Service

CERTIFICATE

No. Z1 10 11 22430 023

Holder of Certificate: **Alu Rehab AS**

Bedriftsveien 23
4353 Klepp Stasjon
NORWAY

Production

57370

Facility(ies):

Certification Mark:



Product:

**Standard wheelchairs
Manual wheelchairs**

Model(s):

Type: Netti 4U Comfort CED/CEDS

Parameters:

Netti 4U comfort CEDS

total length	
incl. footrests:	1110 mm
total length	
without footrests:	840 mm
width:	530 mm – 780 mm
height:	1130 mm

Netti 4U comfort CED

total length	
incl. footrests:	1160 mm
total length	
without footrests:	890 mm
width:	530 mm – 780 mm
height:	1130 mm

Remark: further parameters see enclosure

Tested according to: DIN EN 12183:2009

The product was tested on a voluntary basis and complies with the essential requirements. The certification mark shown above can be affixed on the product. It is not permitted to alter the certification mark in any way. In addition the certification holder must not transfer the certificate to third parties. See also notes overleaf.

Test report no.: 71375630

Date, 2010-11-23

(Torsten Zimmer)

Page 1 of 2





Product Service

Enclosure to certificate no.: Z1 10 11 22430 023

Model identification for the manual wheelchairs
Type: Netti 4U comfort CED/CEDS

	Netti 4U comfort CED(S)	Netti 4U comfort CED
total length incl. footrests:	1110 mm	1160 mm
total length without footrests:	840 mm	890 mm
width:	530 - 780 mm	
height:	1130 mm	
seat width:	350, 400, and 430 mm	430, 450, 500, 550 and 600 mm
seat depth:	375 - 500 mm	425 - 550 mm
seat height:	465 - 500 mm	
backrest height:	480 mm (600 mm with back extension)	
seat plane angle:	- 5° - +20°	
backrest angle:	92° - 137°	
max. load:	130 kg	160 kg
brake system:	knee lever brake (user brake) drum brake (assistant brake)	
back wheels: - size - tyres - adjustable - fixation	12 x 1 3/4", 16 x 1 3/4", 20 x 1", 22 x 1", 24 x 1" and 24 x 1 3/8" solid rubber (accident proof) yes Quick release	
front wheels: - size - tyres - adjustable - fixation	ø 125 x 27, ø 150 x 27 and 175 x 33 mm solid rubber (accident proof) yes Quick release	
frame material:	aluminium	
armrest:	height adjustable and swing able	
legrest:	height adjustable, angle adjustable and removable	
anti tip device:	adjustable	

Department: MHS5 / me
 Date: 2010-11-23