



Product Service

# CERTIFICATE

No. Z1 022430 0056 Rev. 00

**Holder of Certificate:** **Alu Rehab AS**



Bedriftsvegen 23  
4353 Klepp Stasjon  
NORWAY

**Certification Mark:**



**Product:** **Care wheelchairs**

The product was tested on a voluntary basis and complies with the essential requirements. The certification mark shown above can be affixed on the product. It is not permitted to alter the certification mark in any way. In addition, the certification holder must not transfer the certificate to third parties. This certificate is valid until the listed date, unless it is cancelled earlier. All applicable requirements of the testing and certification regulations of TÜV SÜD Group have to be complied. For details see: [www.tuvsud.com/ps-cert](http://www.tuvsud.com/ps-cert)

**Test report no.:** 713150838

**Valid until:** 2025-07-07

**Date,** 2020-07-10

( Torsten Zimmer )



Product Service

# CERTIFICATE

No. Z1 022430 0056 Rev. 00

**Model(s):** Netti Dynamic S

**Brand Name:** Netti by Alu Rehab

## Parameters:

total length incl. footrests:	905 - 1000 mm
total length w/o footrests:	740 - 835 mm
width:	405 - 515 mm
height:	885 mm (w/o headrest)
seat width:	250, 300 and 350 mm
seat depth:	250 - 350 mm
seat height:	440 / 470 mm
backrest height:	375 mm
seat plane angle:	-2° to +35°
backrest angle:	89° to 125°
max. user weight:	75 kg
brake system:	Drum brake
back wheels:	
- size	16 x 1 3/8" and 22 x 1"
- tyres	Solid rubber (puncture proof)
- adjustable	Yes
- fixation	rigid
front wheels:	
- size	Ø 150 x 30 mm
- tyres	Solid rubber (puncture proof)
- adjustable	Yes
- fixation	rigid
frame material:	aluminum/steel
armrest:	height adjustable
legrest:	height adjustable, swingable and removable
anti tip device:	available

**Tested according to:** DIN EN 12183:2014

**Production** 057370

**Facility(ies):**