



Product Service

CERTIFICATE

No. Z1 022430 0057 Rev. 00

Holder of Certificate: **Alu Rehab AS**



Bedriftsvegen 23
4353 Klepp Stasjon
NORWAY

Certification Mark:



Product: **Care wheelchairs**

The product was tested on a voluntary basis and complies with the essential requirements. The certification mark shown above can be affixed on the product. It is not permitted to alter the certification mark in any way. In addition, the certification holder must not transfer the certificate to third parties. This certificate is valid until the listed date, unless it is cancelled earlier. All applicable requirements of the testing and certification regulations of TÜV SÜD Group have to be complied. For details see: www.tuvsud.com/ps-cert

Test report no.: 713150838-1

Valid until: 2025-07-07

Date, 2020-07-10

(Torsten Zimmer)



Product Service

CERTIFICATE

No. Z1 022430 0057 Rev. 00

Model(s): Netti S

Brand Name: Netti by Alu Rehab

Parameters:

total length incl. footrests:	905 - 1000 mm
total length w/o footrests:	740 - 835 mm
width:	405 - 515 mm
height:	885 mm (w/o headrest)
seat width:	250, 300 and 350 mm
seat depth:	250 - 350 mm
seat height:	440 / 470 mm
backrest height:	375 mm
seat plane angle:	-2° to +35°
backrest angle:	89° to 125°
max. user weight:	75 kg
brake system:	Drum brake
back wheels:	
- size	16 x 1 3/8" and 22 x 1"
- tyres	Solid rubber (puncture proof)
- adjustable	Yes
- fixation	rigid
front wheels:	
- size	Ø 150 x 30 mm
- tyres	Solid rubber (puncture proof)
- adjustable	Yes
- fixation	rigid
frame material:	aluminum/steel
armrest:	height adjustable
legrest:	height adjustable, swingable and removable
anti tip device:	available

Tested according to: DIN EN 12183:2014

Production 057370

Facility(ies):