



Combination agreement

The below mentioned ALBER products are in principle appropriate to be combined with the below mentioned ALU REHAB products provided that the combination is made expertly and according to the guidelines of both manufacturers. The guidelines can be asked at the manufacturer.

ALBER product	ALU REHAB products
e-fix E35 (max. 120 kg) / e-fix E36 (max. 160 kg)	Netti II Netti III Netti III EL Netti III HD Netti 4U CED Netti 4U CE Netti 4U CE Plus Netti 4U Base

Please note that mounting has to be effected according to the ALBER mounting instructions by authorized and skilled staff only. Restrictions (seat width, seat height,...) from the ALBER bracket list have to be observed.

ALBER GMBH
Vor dem Weißen Stein 21
72461 Albstadt
Germany

ALU REHAB A/S
Bedriftsveien 23
4353 Klepp Stasjon
Norway

Albstadt, 2018-10-23

Klepp Stasjon,

Hansjörg Reiner
Managing Director

Sigrun Peerstøe-Kotthaus
Product safety responsible